

ABX Minotrol Retic



LOT RX 446

CONTROL

(Exp.) **2024-05-05**
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL			TOLERANCES TOLERANCE	CONTROL			TOLERANCES TOLERANCE	CONTROL			TOLERANCES TOLERANCE
		YUMIZEN		1		YUMIZEN		2		YUMIZEN		3	
		H2500 HELO 2.0				H2500 HELO 2.0				H2500 HELO 2.0			
Méthode Manuelle / Manual													
RET	%	1.6			± 0.8	5.4			± 2.2	10.4			± 3.2
RET	%	1.63			± 0.82	4.71			± 1.18	9.41			± 1.88
RET	10 ⁶ /mm ³	0.057			± 0.029	0.166			± 0.042	0.324			± 0.065
	10 ⁹ /l	57			± 29	166			± 42	324			± 65
CRC	%	1.11			± 0.56	3.20			± 0.96	6.01			± 1.50
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.51			± 0.30	3.52			± 0.30	3.44			± 0.30
IRF	--	0.21			± 0.13	0.26			± 0.16	0.28			± 0.17
RHCC	pg	23.1			± 4.6	16.0			± 3.2	11.2			± 2.8

Ref: TEMP-0382 Rev.45 1300116193 BACK / VERSO