

ABX Minocal



LOT

CX 493

CAL



(Exp.) 2024-05-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN		YUMIZEN						
		H1500 H2500 HELO 1.x	H1500 H2500 HELO 2.0							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.65	8.65							± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.53	4.56							± 0.06
HB HGB	g/dl	13.5	13.5							± 0.2
	g/l	135	135							± 2
	mmol/l	8.38	8.38							± 0.12
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	240	238							± 10
VGM MCV	µm ³ ; fl	83.0	83.9							± 4.0

1300116192
Ref: TEMP-0387 Rev.51