

ABX Minocal



LOT

CX 494

CAL



(Exp.) 2024-06-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio								TOLERANCES TOLERANCE
		YUMIZEN		YUMIZEN						
		H1500 H2500 HELO 1.x	H1500 H2500 HELO 2.0							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.60	8.55							± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52	4.53							± 0.06
HB HGB	g/dl	13.3	13.5							± 0.2
	g/l	133	135							± 2
	mmol/l	8.26	8.38							± 0.12
PLA PLT	10 ³ /mm ³ ;10 ⁹ /l	243	243							± 10
VGM MCV	µm ³ ; fl	83.0	84.5							± 4.0

1300116192
Ref: TEMP-0387 Rev.51