

ABX Minotrol Retic



LOT RX 447
Rev 1

CONTROL

(Exp.) 2024-07-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL			TOLERANCES TOLERANCE	CONTROL			TOLERANCES TOLERANCE	CONTROL			TOLERANCES TOLERANCE
		YUMIZEN		1		YUMIZEN		2		YUMIZEN		3	
		H2500 HELO 1.x				H2500 HELO 1.x				H2500 HELO 1.x			
Méthode Manuelle / Manual													
RET	%	1.4			± 0.7	4.8			± 2.0	10.7			± 3.3
RET	%	1.81			± 0.91	5.22			± 1.31	9.92			± 1.98
RET	10 ⁶ /mm ³	0.062			± 0.031	0.183			± 0.046	0.349			± 0.070
	10 ⁹ /l	62			± 31	183			± 46	349			± 70
MFI MRI *	%	24.4			± 8.0	30.2			± 10.0	31.4			± 12.0
	-	0.244			± 0.080	0.302			± 0.100	0.314			± 0.120
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.43			± 0.30	3.50			± 0.30	3.52			± 0.30
PIC RANGE	-	19			± 8	19			± 8	19			± 8

Ref: TEMP-0382 Rev.45 FRONT / RECTO 1300116193

* For research use only in Canada and USA