

ABX Minotrol Retic



LOT RX 447

CONTROL

(Exp.) **2024-07-05**
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX											
		CONTROL 1			TOLERANCES TOLERANCE	CONTROL 2			TOLERANCES TOLERANCE	CONTROL 3			TOLERANCES TOLERANCE
		YUMIZEN				YUMIZEN				YUMIZEN			
		H2500 HELO 2.0				H2500 HELO 2.0				H2500 HELO 2.0			
Méthode Manuelle / Manual													
RET	%	1.4			± 0.7	4.8			± 2.0	10.7			± 3.3
RET	%	1.81			± 0.91	5.22			± 1.31	9.92			± 1.98
RET	10 ⁶ /mm ³	0.062			± 0.031	0.182			± 0.046	0.349			± 0.070
	10 ⁹ /l	62			± 31	182			± 46	349			± 70
CRC	%	1.19			± 0.60	3.48			± 1.04	6.55			± 1.64
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	3.43			± 0.30	3.49			± 0.30	3.52			± 0.30
IRF	--	0.25			± 0.15	0.31			± 0.19	0.34			± 0.20
RHCC	pg	20.9			± 4.2	15.0			± 3.0	10.7			± 2.7

Ref: TEMP-0382 Rev.45 1300116193 BACK / VERSO