

MINOTROL CRP

Hematology Reference Control

9755 Research Drive
 Irvine, California 92618
 Telephone: (949) 453-0500
 (888) 903-5001
 Fax: (949) 753-0533

Lot No.: **MC402**
 Exp. Date: **2017-04-05** (Y/M/D)

	Parameter/Units	Low Abnormal Lot # MC402-1			Normal Lot # MC402-2			High Abnormal Lot # MC402-3		
		Assay	Range	Limits	Assay	Range	Limits	Assay	Range	Limits
MICROS CRP 200	WBC K/ μ L	2.0	0.4	1.6 - 2.4	7.6	1.0	6.6 - 8.6	20.1	2.2	17.9 - 22.3
	RBC M/ μ L	2.42	0.15	2.27 - 2.57	4.72	0.18	4.54 - 4.90	5.65	0.25	5.40 - 5.90
	HGB g/dL	6.3	0.5	5.8 - 6.8	13.7	0.6	13.1 - 14.3	17.8	0.7	17.1 - 18.5
	HCT %	16.9	2.0	14.9 - 18.9	37.8	2.5	35.3 - 40.3	50.3	3.0	47.3 - 53.3
	MCV fL	70	4	66 - 74	80	4	76 - 84	89	4	85 - 93
	MCH pg	26.0	2.4	23.6 - 28.4	29.0	2.4	26.6 - 31.4	31.5	2.8	28.7 - 34.3
	MCHC g/dL	37.2	3.0	34.2 - 40.2	36.3	3.0	33.3 - 39.3	35.4	3.0	32.4 - 38.4
	RDW %	13.8	4.0	9.8 - 17.8	13.2	4.0	9.2 - 17.2	12.8	4.0	8.8 - 16.8
	PLT K/ μ L	70	20	50 - 90	255	40	215 - 295	510	60	450 - 570
	MPV fL	9.0	2.0	7.0 - 11.0	8.7	2.0	6.7 - 10.7	8.6	2.0	6.6 - 10.6
	Lym %	64.0	8.0	56.0 - 72.0	34.0	6.0	28.0 - 40.0	16.5	6.0	10.5 - 22.5
	Lym K/ μ L	1.3	0.3	1.0 - 1.6	2.6	0.6	2.0 - 3.2	3.3	1.3	2.0 - 4.6
	Mon %	10.0	7.0	3.0 - 17.0	7.0	5.0	2.0 - 12.0	5.5	5.0	0.5 - 10.5
	Mon K/ μ L	0.2	0.2	0.0 - 0.4	0.5	0.4	0.1 - 0.9	1.1	1.0	0.1 - 2.1
	Gran %	26.0	8.0	18.0 - 34.0	59.0	6.0	53.0 - 65.0	78.0	7.0	71.0 - 85.0
	Gran K/ μ L	0.5	0.3	0.2 - 0.8	4.5	0.6	3.9 - 5.1	15.7	1.5	14.2 - 17.2
	CRP mg/dL	0.54	0.30	0.24 - 0.84	2.27	0.70	1.57 - 2.97	4.28	1.00	3.28 - 5.28

For technical assistance, call the Customer Support Center at 1-888-903-5001.