

# MINOTROL CRP

## Hematology Reference Control

9755 Research Drive  
 Irvine, California 92618  
 Telephone: (949) 453-0500  
 (888) 903-5001  
 Fax: (949) 753-0533

Lot No.: **MC403**  
 Exp. Date: **2017-07-05** (Y/M/D)

	Parameter/Units	Low Abnormal Lot # <b>MC403-1</b>			Normal Lot # <b>MC403-2</b>			High Abnormal Lot # <b>MC403-3</b>		
		Assay	Range	Limits	Assay	Range	Limits	Assay	Range	Limits
MICROS CRP 200	WBC K/ $\mu$ L	<b>2.0</b>	0.4	1.6 - 2.4	<b>7.7</b>	1.0	6.7 - 8.7	<b>20.3</b>	2.2	18.1 - 22.5
	RBC M/ $\mu$ L	<b>2.25</b>	0.15	2.10 - 2.40	<b>4.68</b>	0.18	4.50 - 4.86	<b>5.70</b>	0.25	5.45 - 5.95
	HGB g/dL	<b>5.9</b>	0.5	5.4 - 6.4	<b>13.5</b>	0.6	12.9 - 14.1	<b>17.8</b>	0.7	17.1 - 18.5
	HCT %	<b>16.0</b>	2.0	14.0 - 18.0	<b>37.9</b>	2.5	35.4 - 40.4	<b>50.7</b>	3.0	47.7 - 53.7
	MCV fL	<b>71</b>	4	67 - 75	<b>81</b>	4	77 - 85	<b>89</b>	4	85 - 93
	MCH pg	<b>26.2</b>	2.4	23.8 - 28.6	<b>28.8</b>	2.4	26.4 - 31.2	<b>31.2</b>	2.8	28.4 - 34.0
	MCHC g/dL	<b>36.9</b>	3.0	33.9 - 39.9	<b>35.6</b>	3.0	32.6 - 38.6	<b>35.1</b>	3.0	32.1 - 38.1
	RDW %	<b>14.0</b>	4.0	10.0 - 18.0	<b>14.0</b>	4.0	10.0 - 18.0	<b>13.5</b>	4.0	9.5 - 17.5
	PLT K/ $\mu$ L	<b>75</b>	20	55 - 95	<b>265</b>	40	225 - 305	<b>515</b>	60	455 - 575
	MPV fL	<b>9.2</b>	2.0	7.2 - 11.2	<b>8.7</b>	2.0	6.7 - 10.7	<b>8.4</b>	2.0	6.4 - 10.4
	Lym %	<b>62.5</b>	7.0	55.5 - 69.5	<b>34.5</b>	6.0	28.5 - 40.5	<b>16.0</b>	6.0	10.0 - 22.0
	Lym K/ $\mu$ L	<b>1.3</b>	0.3	1.0 - 1.6	<b>2.7</b>	0.6	2.1 - 3.3	<b>3.2</b>	1.3	1.9 - 4.5
	Mon %	<b>9.5</b>	6.0	3.5 - 15.5	<b>8.5</b>	5.0	3.5 - 13.5	<b>6.0</b>	5.0	1.0 - 11.0
	Mon K/ $\mu$ L	<b>0.2</b>	0.2	0.0 - 0.4	<b>0.7</b>	0.4	0.3 - 1.1	<b>1.2</b>	1.0	0.2 - 2.2
	Gran %	<b>28.0</b>	8.0	20.0 - 36.0	<b>57.0</b>	6.0	51.0 - 63.0	<b>78.0</b>	7.0	71.0 - 85.0
	Gran K/ $\mu$ L	<b>0.6</b>	0.3	0.3 - 0.9	<b>4.4</b>	0.6	3.8 - 5.0	<b>15.8</b>	1.5	14.3 - 17.3
	CRP mg/dL	<b>0.76</b>	0.30	0.46 - 1.06	<b>2.49</b>	0.70	1.79 - 3.19	<b>3.94</b>	1.00	2.94 - 4.94

For technical assistance, call the Customer Support Center at 1-888-903-5001.