

MINOTROL CRP

Hematology Reference Control

9755 Research Drive
 Irvine, California 92618
 Telephone: (949) 453-0500
 (888) 903-5001
 Fax: (949) 753-0533

Lot No.: **MC404**
 Exp. Date: **2017-10-05** (Y/M/D)

| | Parameter/Units | Low Abnormal Lot # MC404-1 | | | Normal Lot # MC404-2 | | | High Abnormal Lot # MC404-3 | | |
|----------------|-----------------|-----------------------------------|-------|-------------|-----------------------------|-------|-------------|------------------------------------|-------|-------------|
| | | Assay | Range | Limits | Assay | Range | Limits | Assay | Range | Limits |
| MICROS CRP 200 | WBC K/ μ L | 2.0 | 0.4 | 1.6 - 2.4 | 7.4 | 1.0 | 6.4 - 8.4 | 19.8 | 2.2 | 17.6 - 22.0 |
| | RBC M/ μ L | 2.43 | 0.15 | 2.28 - 2.58 | 4.77 | 0.18 | 4.59 - 4.95 | 5.76 | 0.25 | 5.51 - 6.01 |
| | HGB g/dL | 6.2 | 0.5 | 5.7 - 6.7 | 13.4 | 0.6 | 12.8 - 14.0 | 17.9 | 0.7 | 17.2 - 18.6 |
| | HCT % | 17.0 | 2.0 | 15.0 - 19.0 | 38.2 | 2.5 | 35.7 - 40.7 | 51.3 | 3.0 | 48.3 - 54.3 |
| | MCV fL | 70 | 4 | 66 - 74 | 80 | 4 | 76 - 84 | 89 | 4 | 85 - 93 |
| | MCH pg | 25.5 | 2.4 | 23.1 - 27.9 | 28.1 | 2.4 | 25.7 - 30.5 | 31.1 | 2.8 | 28.3 - 33.9 |
| | MCHC g/dL | 36.4 | 3.0 | 33.4 - 39.4 | 35.1 | 3.0 | 32.1 - 38.1 | 34.9 | 3.0 | 31.9 - 37.9 |
| | RDW % | 14.0 | 4.0 | 10.0 - 18.0 | 14.5 | 4.0 | 10.5 - 18.5 | 13.5 | 4.0 | 9.5 - 17.5 |
| | PLT K/ μ L | 75 | 20 | 55 - 95 | 250 | 40 | 210 - 290 | 485 | 60 | 425 - 545 |
| | MPV fL | 9.0 | 2.0 | 7.0 - 11.0 | 8.3 | 2.0 | 6.3 - 10.3 | 7.6 | 2.0 | 5.6 - 9.6 |
| | Lym % | 65.5 | 7.0 | 58.5 - 72.5 | 33.5 | 6.0 | 27.5 - 39.5 | 15.5 | 6.0 | 9.5 - 21.5 |
| | Lym K/ μ L | 1.3 | 0.2 | 1.1 - 1.5 | 2.5 | 0.5 | 2.0 - 3.0 | 3.1 | 1.2 | 1.9 - 4.3 |
| | Mon % | 8.5 | 6.0 | 2.5 - 14.5 | 8.5 | 5.0 | 3.5 - 13.5 | 5.5 | 4.5 | 1.0 - 10.0 |
| | Mon K/ μ L | 0.2 | 0.2 | 0.0 - 0.4 | 0.6 | 0.4 | 0.2 - 1.0 | 1.1 | 0.9 | 0.2 - 2.0 |
| | Gran % | 26.0 | 8.0 | 18.0 - 34.0 | 58.0 | 6.0 | 52.0 - 64.0 | 79.0 | 7.0 | 72.0 - 86.0 |
| | Gran K/ μ L | 0.5 | 0.2 | 0.3 - 0.7 | 4.3 | 0.5 | 3.8 - 4.8 | 15.6 | 1.4 | 14.2 - 17.0 |
| | CRP mg/dL | 0.54 | 0.30 | 0.24 - 0.84 | 2.32 | 0.70 | 1.62 - 3.02 | 4.78 | 1.00 | 3.78 - 5.78 |

For technical assistance, call the Customer Support Center at 1-888-903-5001.