



**HORIBA ABX SAS**  
 Parc Euromédecine - Rue du Caducée  
 B.P. 7290  
 34184 Montpellier Cedex 4  
 FRANCE  
 www.horiba.com



**LOT** 1: 1706101



1 7 0 6 1 0 1

**REF** : 1300026187



2019/04/30

1300026587



1300026587

(FR) Constituant (EN) Constituent (DE) Bestandteil (ES) Componente (IT) Componente (PT) Componente (DA) Bestanddel (SV) Konstituent (EL) συστατικό (PL) Składnik (NO) Bestanddel	ABX Pentra 400	Pentra C400	Pentra C200	Unité	Valeur titrée	Intervalle de confiance	
				Unit Einheit Unidad Unità Unidade Enhed Enhet Μονάδα Jednostka Enhet	Assigned value Deklariertes Wert Valor asignado Valori assegnati Valor determinado Tildelt værdi Tildelat värde Αναγραφόμενη τιμή Wartość przypisana Tildelt verd	Confidence range Vertrauensintervalle Intervalos de confianza Livelli di fiducia Intervalo de confiança Konfidensområde Konfidensintervall Εύρος εμπιστοσύνης Zakres ufności Konfidensspekter	1 SD
C				U	$\oplus$	Confidence range	1 SD
CK (IFCC/DGKC) (37°C)	●			U/L μkat/L	<b>167</b> <b>2.79</b>	149 - 185 2.49 - 3.09	6 0.10
		●		U/L μkat/L	<b>167</b> <b>2.79</b>	149 - 185 2.49 - 3.09	6 0.10
			●	U/L μkat/L	<b>164</b> <b>2.74</b>	149 - 179 2.50 - 2.98	5 0.08
CK-MB (Immunology method) (37°C)	●			U/L μkat/L	<b>34.8</b> <b>0.58</b>	29.1 - 40.5 0.49 - 0.67	1.9 0.03
		●		U/L μkat/L	<b>34.8</b> <b>0.58</b>	29.1 - 40.5 0.49 - 0.67	1.9 0.03
			●	U/L μkat/L	<b>41.9</b> <b>0.70</b>	35.6 - 48.2 0.58 - 0.82	2.1 0.04