



Yumizen G PT 5 for Yumizen G Line



A 4 9 7 0 1 3 0 3 6 3 3 8 0 1 1 9 A



A 5 1 0 0 1 3 5 5 0 0 2 0 1 3 3 0 A



A 6 2 6 8 2 5 0 3 3 2 1 3 5 1 1 3 A



A 7 1 4 5 A

Reference No.:		1300036338		
Lot No.:		970130		
Expiry date:		2019-01-31		
Coagulometers		ISI	MASTER CURVE	
Optical	Yumizen G Line		%	sec
		100	13,5	
		50	20,1	
		33	26,8	
		25	33,2	
		1,13	MNPT	13,5
Important notice				
<ul style="list-style-type: none"> - The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT. - The efficiency of the reconstitution of the reagent may be approved by performing it at 37°C. 				