



Yumizen G PT 5 for Yumizen G Line



A 4 9 8 0 5 1 4 3 6 3 3 8 0 5 2 0 A



A 5 1 0 0 1 3 9 5 0 0 2 0 0 3 3 0 A



A 6 2 6 3 2 5 0 3 2 3 1 3 9 1 1 4 A



A 7 1 5 1 A

Reference No.:		1300036338		
Lot No.:		980514		
Expiry date:		2020-05-31		
Coagulometers		ISI	MASTER CURVE	
Optical	Yumizen G Line		1,14	%
		100		13,9
		50		20,0
		33		26,3
		25		32,3
			MNPT	13,9
Important notice				
<ul style="list-style-type: none"> - The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT. - The efficiency of the reconstitution of the reagent may be approved by performing it at 37°C. 				