



Yumizen G PT Liq 4 for Yumizen G Line



A 4 9 7 0 6 1 3 3 6 3 7 3 0 8 1 8 A



A 5 1 0 0 1 2 7 5 0 0 1 8 2 3 3 0 A



A 6 2 3 9 2 5 0 2 9 3 1 2 7 1 1 7 A



A 7 1 7 4 A

Reference No.:		1300036373		
Lot No.:		970613		
Expiry date:		2018-08-31		
Coagulometers		ISI	MASTER CURVE	
Optical	Yumizen G Line		%	sec
		100	12,7	
		50	18,2	
		33	23,9	
		25	29,3	
		1,17	MNPT 12,7	
Important notice				
<p>The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.</p>				