



# ***Yumizen G FIB 2*** ***for Yumizen G Line***



|                       |                       |                     |             |
|-----------------------|-----------------------|---------------------|-------------|
| <b>Reference No.:</b> |                       | <b>1300036383</b>   |             |
| <b>Lot No.:</b>       |                       | <b>960805</b>       |             |
| <b>Expiry date:</b>   |                       | <b>2018-08-31</b>   |             |
| <b>Coagulometers</b>  |                       | <b>MASTER CURVE</b> |             |
| <b>Optical</b>        | <b>Yumizen G Line</b> | <b>g/L</b>          | <b>sec</b>  |
|                       |                       | <b>4</b>            | <b>7,2</b>  |
|                       |                       | <b>3</b>            | <b>9,2</b>  |
|                       |                       | <b>2</b>            | <b>13,1</b> |
|                       |                       | <b>1</b>            | <b>23,6</b> |
| <b>Dilution</b>       |                       | <b>10x</b>          |             |