



Yumizen G CTRL I (level 1) for Yumizen G Line



Reference No.:		1300036412	
Lot No:		970614_REV2	
Expiry date:		2019-06-30	
Optical	Yumizen G Line		
Kit	Parameter	Unit	Range
PT	Prothrombin Time	sec	10,4-14,1
	Prothrombin %	%	83-139
	Prothrombin INR	INR	0,76-1,14
PT Liq	Prothrombin Time	sec	10,9-14,7
	Prothrombin %	%	86-144
	Prothrombin INR	INR	0,78-1,18
PT R	Prothrombin Time	sec	8,4-11,4
	Prothrombin %	%	93-155
	Prothrombin INR	INR	0,76-1,14
APTT	Partial Thromboplastin Time	sec	32,1-43,5
APTT Liq	Partial Thromboplastin Time	sec	28,9-39,1
FIB Cl	Fibrinogen	g/L	1,9-2,9
TT	Thrombin Time	sec	20,1-30,2



Yumizen G CTRL II (level 2) for Yumizen G Line



D 4 9 7 0 6 1 4 5 6 4 1 2 0 6 1 9 D



D 5 1 5 8 2 3 7 0 3 3 0 5 5 1 3 4 D



D 6 2 0 1 1 5 4 2 4 4 0 3 6 0 6 0 D



D 7 1 2 7 1 9 1 1 1 6 1 7 4 0 4 1 D



D 8 0 6 9 1 1 5 1 7 2 5 5 3 7 5 3 D



D 9 4 9 6 6 7 0 0 7 0 1 6 0 0 0 0 D



D 0 9 9 9 3 6 7 D

Reference No.:		1300036412	
Lot No:		970614_REV2	
Expiry date:		2019-06-30	
Optical	Yumizen G Line		
Kit	Parameter	Unit	Range
PT	Prothrombin Time	sec	15,8-23,7
	Prothrombin %	%	33-55
	Prothrombin INR	INR	1,34-2,01
PT Liq	Prothrombin Time	sec	15,4-24,4
	Prothrombin %	%	36-60
	Prothrombin INR	INR	1,27-1,91
PT R	Prothrombin Time	sec	11,6-17,4
	Prothrombin %	%	41-69
	Prothrombin INR	INR	1,15-1,72
APTT	Partial Thromboplastin Time	sec	55,3-75,3
APTT Liq	Partial Thromboplastin Time	sec	49,6-67,0
FIB Cl	Fibrinogen	g/L	0,7-1,6
TT	Thrombin Time	sec	NA