



# ***Yumizen G CTRL I (level 1) for Yumizen G Line***



D4980213464120220D



D5100140088146077D



D6115100140083139D



D7076114083112095D



D8176068102310419D



D9290393200320212D



D0318314D

Reference No.:		1300036412	
Lot No:		980213_REV2	
Expiry date:		2020-02-29	
Optical	Yumizen G Line		
Kit	Parameter	Unit	Range
PT	Prothrombin Time	sec	10,0-14,0
	Prothrombin %	%	88-146
	Prothrombin INR	INR	0,77-1,15
PT Liq	Prothrombin Time	sec	10,0-14,0
	Prothrombin %	%	83-139
	Prothrombin INR	INR	0,76-1,14
PT R	Prothrombin Time	sec	8,3-11,2
	Prothrombin %	%	95-176
	Prothrombin INR	INR	0,68-1,02
APTT	Partial Thromboplastin Time	sec	31,0-41,9
APTT Liq	Partial Thromboplastin Time	sec	29,0-39,3
FIB Cl	Fibrinogen	g/L	2,0-3,2
TT	Thrombin Time	sec	21,2-31,8
AT	Antithrombin III	%	84-114



# Yumizen G CTRL II (level 2) for Yumizen G Line



D 4 9 8 0 2 1 3 5 6 4 1 2 0 2 2 0 D



D 5 1 4 5 2 1 8 0 3 9 0 6 6 1 2 0 D



D 6 1 8 0 1 4 3 2 3 3 0 3 5 0 5 8 D



D 7 1 3 0 1 9 6 1 1 3 1 6 9 0 4 6 D



D 8 0 8 5 0 9 9 1 4 9 5 0 3 7 0 3 D



D 9 4 9 1 6 6 4 0 8 0 1 9 0 0 0 0 D



D 0 9 9 9 3 6 7 D

Reference No.:		1300036412	
Lot No:		980213_REV2	
Expiry date:		2020-02-29	
Optical	Yumizen G Line		
Kit	Parameter	Unit	Range
PT	Prothrombin Time	sec	14,5-21,8
	Prothrombin %	%	39-66
	Prothrombin INR	INR	1,20-1,80
PT Liq	Prothrombin Time	sec	14,3-23,3
	Prothrombin %	%	35-58
	Prothrombin INR	INR	1,30-1,96
PT R	Prothrombin Time	sec	11,3-16,9
	Prothrombin %	%	46-85
	Prothrombin INR	INR	0,99-1,49
APTT	Partial Thromboplastin Time	sec	50,3-70,3
APTT Liq	Partial Thromboplastin Time	sec	49,1-66,4
FIB Cl	Fibrinogen	g/L	0,8-1,9
TT	Thrombin Time	sec	NA
AT	Antithrombin III	%	36-61