



# Yumizen G PT Liq 4 for Yumizen G Line



<b>Reference No.:</b>		<b>1300036373</b>			
<b>Lot No.:</b>		<b>980813</b>			
<b>Expiry date:</b>		<b>2020-02-29</b>			
<b>Coagulometers</b>		<b>ISI</b>	<b>MASTER CURVE</b>		
			<b>%</b>	<b>sec</b>	
<b>Optical</b>	<b>Yumizen G Line</b>	<b>1,18</b>	<b>100</b>	<b>12,6</b>	
			<b>50</b>	<b>18,0</b>	
			<b>33</b>	<b>23,6</b>	
			<b>25</b>	<b>28,9</b>	
			<b>MNPT</b>	<b>12,6</b>	
<b>Important notice</b>					
The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.					