



Yumizen G DDi 2 for YUMIZEN G1500



Reference No.:		1300036391	
Lot No.:		981114	
Expiry date:		2019-11-30	
Coagulometers		MASTER CURVE	
		µg(FEU)/mL	ΔOD
Optical	Yumizen G1500/G800	4,00	0,3840
		3,00	0,2860
		2,00	0,1420
		1,00	0,0575
		0,50	0,0285
		0,25	0,0162
Optical	Yumizen G400/G200	4,00	0,2830
		3,00	0,2300
		2,00	0,1540
		1,00	0,0680
		0,50	0,0340
		0,25	0,0190