



HORIBA ABX SAS
 Parc Euromédecine - Rue du Caducée
 B.P. 7290
 34184 Montpellier Cedex 4
 FRANCE
 www.horiba.com



LOT 1: 1902401



1 9 0 2 4 0 1

REF : 1300026187



2021/07/31

SAP 1300056305



1300056305

(FR) Constituant (EN) Constituent (DE) Bestandteil (ES) Componente (IT) Componente (PT) Componente (DA) Bestanddel (SV) Konstituent (EL) συστατικό (PL) Składnik (NO) Bestanddel	ABX Pentra 400	Pentra C400	Pentra C200	Unité Unit Einheit Unidad Unità Unidade Enhed Enhet Μονάδα Jednostka Enhet	Valeur titrée Assigned value Deklariertes Wert Valor asignado Valori assegnati Valor determinado Tildelt værdi Tildelat värde Αναγραφόμενη τιμή Wartość przypisana Tildelt verd	Intervalle de confiance Confidence range Vertrauensintervalle Intervalos de confianza Livelli di fiducia Intervalo de confiança Konfidensområde Konfidensintervall Εύρος εμπιστοσύνης Zakres ufności Konfidensspekter	
C				U	\oplus	Confidence range*	1 SD
CK (IFCC/DGKC) (37°C)	●			U/L μkat/L	150 2.51	134 - 166 2.25 - 2.77	8 0.13
		●		U/L μkat/L	150 2.51	134 - 166 2.25 - 2.77	8 0.13
			●	U/L μkat/L	153 2.56	137 - 169 2.17 - 2.95	8 0.13
CK-MB (Immunology method) (37°C)	●			U/L μkat/L	33.1 0.550	28.7 - 37.5 0.470 - 0.630	2.2 0.040
		●		U/L μkat/L	33.1 0.550	28.7 - 37.5 0.470 - 0.630	2.2 0.040
			●	U/L μkat/L	42.6 0.711	36.2 - 49.0 0.552 - 0.870	3.2 0.053

* Confidence range: +/- 2 SD