



HORIBA ABX SAS
 BP 7290
 34184 Montpellier
 cedex 4 - France
 www.horiba-
 abx.com



A03A00013A

Yumizen G FIB 5 for Yumizen G Line



D4981009363841020D



D5400062300080200D



D6112100203000000D



D7111D

Reference No.:		1300036384	
Lot No.:		981009	
Expiry date:		2020-10-31	
Coagulometers		MASTER CURVE	
		g/L	sec
Optical	Yumizen G Line	4	6,2
		3	8,0
		2	11,2
		1	20,3
		Dilution	10x