



Yumizen G CTRL I (level 1) for Yumizen G Line



D 4 9 8 0 8 0 7 4 6 4 1 2 0 8 2 0 D



D 5 1 1 6 1 5 7 0 7 4 1 2 3 0 8 1 D



D 6 1 2 2 1 1 5 1 5 5 0 6 9 1 1 5 D



D 7 0 8 4 1 2 6 0 8 8 1 1 9 0 8 7 D



D 8 1 4 6 0 7 5 1 1 2 3 4 0 4 6 0 D



D 9 2 9 2 4 2 2 1 8 0 2 7 0 2 3 5 D



D 0 3 5 2 3 5 6 D

Reference No.:		1300036412	
Lot No:		980807_REV3	
Expiry date:		2020-08-31	
Optical	Yumizen G Line		
Kit	Parameter	Unit	Range
PT	Prothrombin Time	sec	11,6-15,7
	Prothrombin %	%	74-123
	Prothrombin INR	INR	0,81-1,22
PT Liq	Prothrombin Time	sec	11,5-15,5
	Prothrombin %	%	69-115
	Prothrombin INR	INR	0,84-1,26
PT R	Prothrombin Time	sec	8,8-11,9
	Prothrombin %	%	87-146
	Prothrombin INR	INR	0,75-1,12
APTT	Partial Thromboplastin Time	sec	34,0-46,0
➔ APTT Liq	Partial Thromboplastin Time	sec	29,2-42,2
FIB Cl	Fibrinogen	g/L	1,8-2,7
TT	Thrombin Time	sec	23,5-35,2
AT	Antithrombin III	%	81-109

→ modification since the previous version



Yumizen G CTRL II (level 2) for Yumizen G Line



D 4 9 8 0 8 0 7 5 6 4 1 2 0 8 2 0 D



D 5 1 6 9 2 5 4 0 3 4 0 5 6 1 3 6 D



D 6 2 0 5 1 6 9 2 5 3 0 2 8 0 4 7 D



D 7 1 4 5 2 1 8 1 2 0 1 8 0 0 4 4 D



D 8 0 7 4 1 1 3 1 6 9 5 5 7 7 5 3 D



D 9 4 8 5 6 9 7 0 7 0 1 7 0 0 0 0 D



D 0 9 9 9 3 9 4 D

Reference No.:		1300036412	
Lot No:		980807_REV3	
Expiry date:		2020-08-31	
Optical	Yumizen G Line		
Kit	Parameter	Unit	Range
PT	Prothrombin Time	sec	16,9-25,4
	Prothrombin %	%	34-56
	Prothrombin INR	INR	1,36-2,05
PT Liq	Prothrombin Time	sec	16,9-25,3
	Prothrombin %	%	28-47
	Prothrombin INR	INR	1,45-2,18
PT R	Prothrombin Time	sec	12,0-18,0
	Prothrombin %	%	44-74
	Prothrombin INR	INR	1,13-1,69
APTT	Partial Thromboplastin Time	sec	55,7-75,3
➔ APTT Liq	Partial Thromboplastin Time	sec	48,5-69,7
FIB Cl	Fibrinogen	g/L	0,7-1,7
TT	Thrombin Time	sec	NA
AT	Antithrombin III	%	34-56

→ modification since the previous version