



Yumizen G PT Liq 4 for Yumizen G Line



Reference No.:		1300036373		
Lot No.:		990724		
Expiry date:		2021-01-31		
Coagulometers		ISI	MASTER CURVE	
			%	sec
Optical	Yumizen G Line	1,29	100	13,2
			50	19,2
			33	25,3
			25	31,1
			MNPT	13,2
Important notice				
- The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.				