



Yumizen G PT Reco 10 for Yumizen G Line



A 4 9 9 0 8 0 5 3 6 3 7 6 0 8 2 1 A



A 5 1 0 0 1 1 2 5 0 0 1 5 0 3 3 0 A



A 6 1 9 0 2 5 0 2 2 7 1 1 2 0 9 4 A



A 7 1 5 6 A

Reference No.:		1300036376	
Lot No.:		990805	
Expiry date:		2021-08-31	
Coagulometers		ISI	MASTER CURVE
Optical	Yumizen G Line		%
		100	11,2
		50	15,0
		33	19,0
		25	22,7
		0,94	MNPT
			11,2
Important notice			
The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.			