



Yumizen G PT Reco 10 for Yumizen G Line


 A4991107363761121A


 A5100096500129330A


 A6163250196096097A


 A7182A

Reference No.:		1300036376	
Lot No.:		991107	
Expiry date:		2021-11-30	
Coagulometers		ISI	MASTER CURVE
Optical	Yumizen G Line		%
		100	9,6
		50	12,9
		33	16,3
		25	19,6
	0,97	MNPT	9,6
Important notice			
The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.			