



Yumizen G DDi 2

for YUMIZEN

G1500

G400

G800

for Yumizen G1500



D4991201363921220D



D5400283030021902D



D6001330100057005D



D7002750250150158D

for Yumizen G800



D4991201363921220D



D5400266530020452D



D6001245100053005D



D7002550250135165D

Reference No.:		1300036391 1300036392	
Lot No.:		991201	
Expiry date:		2020-12-31	
Coagulometers		MASTER CURVE	
		µg(FEU)/mL	ΔOD
Optical I.	Yumizen G1500	4,00	0,2830
		3,00	0,2190
		2,00	0,1330
		1,00	0,0570
		0,50	0,0275
		0,25	0,0150
		Slope value	0,06
Optical II.	Yumizen G400	4,00	0,1148
		3,00	0,0878
		2,00	0,0565
		1,00	0,0263
		0,50	0,0125
		0,25	0,0063
Optical III.	Yumizen G800	4,00	0,2665
		3,00	0,2045
		2,00	0,1245
		1,00	0,0530
		0,50	0,0255
		0,25	0,0135
Slope value	0,06		