

**Yumizen G CTRL DDi I
 for YUMIZEN
 G1500/G1550**



Reference No.:		1300036414
Lot No.:		900224
Expiry date:		2022-02-28
Optical	Yumizen G1500/G1550	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	0,24 – 0,74

**Yumizen G CTRL DDi I
 for YUMIZEN G800**



Reference No.:		1300036414
Lot No.:		900224
Expiry date:		2022-02-28
Optical	Yumizen G800	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	0,27 – 0,77

**Yumizen G CTRL DDi I
 for YUMIZEN
 G400 DDi/G200**

Reference No.:		1300036414
Lot No.:		900224
Expiry date:		2022-02-28
Optical	Yumizen G400 DDi/G200	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	0,23 – 0,73

**Yumizen G CTRL DDi II
 for YUMIZEN
 G1500/G1550**



Reference No.:		1300036414
Lot No.:		900224
Expiry date:		2022-02-28
Optical	Yumizen G1500/G1550	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	1,52 – 2,54

**Yumizen G CTRL DDi II
 for YUMIZEN G800**



Reference No.:		1300036414
Lot No.:		900224
Expiry date:		2022-02-28
Optical	Yumizen G800	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	1,63 – 2,72

**Yumizen G CTRL DDi II
 for YUMIZEN
 G400 DDi/G200**

Reference No.:		1300036414
Lot No.:		900224
Expiry date:		2022-02-28
Optical	Yumizen G400 DDi/G200	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	1,75 – 2,91