



HORIBA ABX SAS
 Parc Euromédecine - Rue du Caducée
 B.P. 7290
 34184 Montpellier Cedex 4
 FRANCE
 www.horiba.com



LOT 1: 2002701 (rev. 1)

REF : 1300026187



2 0 0 2 7 0 1



2023/02/12

SAP 1300056305



1300056305

(FR) Constituant (EN) Constituent (DE) Bestandteil (ES) Componente (IT) Componente (PT) Componente (DA) Bestanddel (SV) Konstituent (EL) συστατικό (PL) Składnik (NO) Bestanddel	ABX Pentra 400	Pentra C400	Pentra C200	Unité	Valeur titrée	Intervalle de confiance	
				Unit	Assigned value	Confidence range	1 SD
				Einheit	Deklariertes Wert	Vertrauensintervalle	
				Unidad	Valor asignado	Intervalos de confianza	
				Unità	Valori assegnati	Livelli di fiducia	
				Unidade	Valor determinado	Intervalo de confiança	
				Enhed	Tildelt værdi	Konfidensområde	
				Enhet	Tildelat värde	Konfidensintervall	
				Μονάδα	Αναγραφόμενη τιμή	Εύρος εμπιστοσύνης	
				Jednostka	Wartość przypisana	Zakres ufności	
				Enhet	Tildelt verd	Konfidensspekter	
C				U	\oplus	Confidence range*	1 SD
CK (IFCC/DGKC) (37°C)	•			U/L	153	137 - 169	8
		•		µkat/L	2.56	2.30 - 2.82	0.13
			•	U/L	153	137 - 169	8
			•	µkat/L	2.56	2.30 - 2.82	0.13
CK-MB (Immunology method) (37°C)	•			U/L	154	138 - 170	8
		•		µkat/L	2.57	2.31 - 2.83	0.13
			•	U/L	30.1	26.1 - 34.1	2.0
	•			µkat/L	0.503	0.437 - 0.569	0.033
		•		U/L	30.1	26.1 - 34.1	2.0
			•	µkat/L	0.503	0.437 - 0.569	0.033
			•	U/L	38.8	33.0 - 44.6	2.9
			•	µkat/L	0.648	0.552 - 0.744	0.048

* Confidence range: +/- 2 SD