

Yumizen G DDi 2 for YUMIZEN G1500/G1550



Reference No.:		1300036391	
Lot No.:		991241 (rev. 1)	
Expiry date:		2020-12-31	
Coagulometers		MASTER CURVE	
Optical	Yumizen G1500/G1550	µg(FEU)/mL	ΔOD
		4,00	0,2480
		3,00	0,1870
		2,00	0,1150
		1,00	0,0510
		0,50	0,0245
		0,25	0,0125
		Slope value	0,06

Yumizen G DDi 2 for YUMIZEN G800



Reference No.:		1300036391	
Lot No.:		991241 (rev. 1)	
Expiry date:		2020-12-31	
Coagulometers		MASTER CURVE	
Optical	Yumizen G800	µg(FEU)/mL	ΔOD
		4,00	0,2305
		3,00	0,1700
		2,00	0,1030
		1,00	0,0455
		0,50	0,0215
		0,25	0,0110
		Slope value	0,06



Yumizen G DDi 2 for YUMIZEN G400 DDi/G200

Reference No.:		1300036391	
Lot No.:		991241 (rev. 1)	
Expiry date:		2020-12-31	
Coagulometers		MASTER CURVE	
		µg(FEU)/mL	ΔOD
Optical	Yumizen G400 DDi/G200	4,00	0,0940
		3,00	0,0710
		2,00	0,0484
		1,00	0,0250
		0,50	0,0128
		0,25	0,0065
		Slope value	N/A
Important notice			
<p>For Yumizen G400 DDi and Yumizen G200, make sure that the end time is set to 180 seconds in the D-Dimer setup and the curve type Lin-Lin p-p is selected in Calibration > D-Dimer > ug FEU/ml. For more information, please contact your local HORIBA Medical representative.</p>			