



Yumizen G DDi 2 for YUMIZEN G1500/G1550



D 4991244363921220D



D 5400275530021452D



D 6001330100056505D



D 7002650250135177D

Reference No.:		1300036391	
Lot No.:		991244 (rev. 1)	
Expiry date:		2020-12-31	
Coagulometers		MASTER CURVE	
Optical	Yumizen G1500/G1550	µg(FEU)/mL	ΔOD
		4,00	0,2755
		3,00	0,2145
		2,00	0,1330
		1,00	0,0565
		0,50	0,0265
		0,25	0,0135
Slope value		0,06	

Yumizen G DDi 2 for YUMIZEN G800



D 4991244363921220D



D 5400262030019902D



D 6001190100050505D



D 7002400250125165D

Reference No.:		1300036391	
Lot No.:		991244 (rev. 1)	
Expiry date:		2020-12-31	
Coagulometers		MASTER CURVE	
Optical	Yumizen G800	µg(FEU)/mL	ΔOD
		4,00	0,2620
		3,00	0,1990
		2,00	0,1190
		1,00	0,0505
		0,50	0,0240
		0,25	0,0125
Slope value		0,06	

Yumizen G DDi 2 ***for YUMIZEN*** ***G400 DDi/G200***

Reference No.:		1300036391	
Lot No.:		991244 (rev. 1)	
Expiry date:		2020-12-31	
Coagulometers		MASTER CURVE	
		µg(FEU)/mL	ΔOD
Optical	Yumizen G400 DDi/G200	4,00	0,1235
		3,00	0,0925
		2,00	0,0580
		1,00	0,0266
		0,50	0,0131
		0,25	0,0070
		Slope value	N/A
Important notice			
<p>For Yumizen G400 DDi and Yumizen G200, make sure that the end time is set to 180 seconds in the D-Dimer setup and the curve type Lin-Lin p-p is selected in Calibration > D-Dimer > ug FEU/ml. For more information, please contact your local HORIBA Medical representative.</p>			