



Yumizen G CTRL I (level 1) for Yumizen G Line



D 4 9 9 0 2 1 1 4 6 4 1 2 0 2 2 1 D



D 5 1 0 1 1 3 6 0 8 5 1 4 1 0 7 5 D



D 6 1 1 3 0 9 8 1 3 3 1 0 4 1 7 3 D



D 7 0 7 0 1 0 5 0 7 6 1 0 3 1 1 4 D



D 8 1 9 0 0 6 8 1 0 1 3 2 8 4 4 4 D



D 9 3 0 3 4 1 0 2 0 0 3 0 0 2 4 9 D



D 0 3 7 4 3 0 3 D

Reference No.:		1300036412	
Lot No:		990211 (rev. 2)	
Expiry date:		2021-02-28	
Optical	Yumizen G Line		
Kit	Parameter	Unit	Range
PT	Prothrombin Time	sec	10,1-13,6
	Prothrombin %	%	85-141
	Prothrombin INR	INR	0,75-1,13
PT Liq	Prothrombin Time	sec	9,8-13,3
	Prothrombin %	%	104-173
	Prothrombin INR	INR	0,70-1,05
PT R	Prothrombin Time	sec	7,6-10,3
	Prothrombin %	%	114-190
	Prothrombin INR	INR	0,68-1,01
APTT	Partial Thromboplastin Time	sec	32,8-44,4
APTT Liq (**)	Partial Thromboplastin Time	sec	30,3-41,0 (**)
FIB Cl	Fibrinogen	g/L	2,0-3,0
TT	Thrombin Time	sec	24,9-37,4
AT	Antithrombin	%	85-115
Important notice			
The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.			
The % and INR dimension are linked to this MNPT. Then, the ranges of the controls need to be adjusted.			

(**):Warning! This range has to be used only with APTTLiq reagents Lots n°: 900412 and higher (90xxxx where xxxx>0412).
For other lots, please refer to the previous revision of values sheet.



Yumizen G CTRL II (level 2) for Yumizen G Line



D 4 9 9 0 2 1 1 5 6 4 1 2 0 2 2 1 D



D 5 1 5 3 2 3 0 0 3 5 0 5 8 1 3 1 D



D 6 1 9 7 1 5 4 2 3 1 0 3 1 0 5 2 D



D 7 1 3 6 2 0 4 1 0 3 1 5 5 0 5 1 D



D 8 0 8 5 0 9 8 1 4 7 5 3 3 7 2 1 D



D 9 4 8 4 6 5 4 0 8 0 1 9 0 0 0 0 D



D 0 9 9 9 3 4 4 D

Reference No.:		1300036412	
Lot No:		990211 (rev. 2)	
Expiry date:		2021-02-28	
Optical	Yumizen G Line		
Kit	Parameter	Unit	Range
PT	Prothrombin Time	sec	15,3-23,0
	Prothrombin %	%	35-58
	Prothrombin INR	INR	1,31-1,97
PT Liq	Prothrombin Time	sec	15,4-23,1
	Prothrombin %	%	31-52
	Prothrombin INR	INR	1,36-2,04
PT R	Prothrombin Time	sec	10,3-15,5
	Prothrombin %	%	51-85
	Prothrombin INR	INR	0,98-1,47
APTT	Partial Thromboplastin Time	sec	53,3-72,1
APTT Liq (**)	Partial Thromboplastin Time	sec	48,4-65,4 (**)
FIB Cl	Fibrinogen	g/L	0,8-1,9
TT	Thrombin Time	sec	NA
AT	Antithrombin	%	37-62
Important notice			
The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.			
The % and INR dimension are linked to this MNPT. Then, the ranges of the controls need to be adjusted.			

(**): Warning! This range has to be used only with APTTLiq reagents Lots n°: 900412 and higher (90xxxx where xxxx>0412).
For other lots, please refer to the previous revision of values sheet.