



Yumizen G DDi 2 for YUMIZEN G1500/G1550



D 4900601363910621D



D 5400274530021302D



D 6001290100053505D



D 7002550250130156D

Reference No.:		1300036391	
Lot No.:		900601 (rev. 1)	
Expiry date:		2021-06-30	
Coagulometers		MASTER CURVE	
Optical	Yumizen G1500/G1550	$\mu\text{g}(\text{FEU})/\text{mL}$	ΔOD
		4,00	0,2745
		3,00	0,2130
		2,00	0,1290
		1,00	0,0535
		0,50	0,0255
		0,25	0,0130
Slope value		0,06	

Yumizen G DDi 2 for YUMIZEN G800



D 4900601363910621D



D 5400256530019652D



D 6001160100048005D



D 7002300250120158D

Reference No.:		1300036391	
Lot No.:		900601 (rev. 1)	
Expiry date:		2021-06-30	
Coagulometers		MASTER CURVE	
Optical	Yumizen G800	$\mu\text{g}(\text{FEU})/\text{mL}$	ΔOD
		4,00	0,2565
		3,00	0,1965
		2,00	0,1160
		1,00	0,0480
		0,50	0,0230
		0,25	0,0120
Slope value		0,06	



Yumizen G DDi 2 for YUMIZEN G400 DDi/G200

Reference No.:		1300036391	
Lot No.:		900601 (rev. 1)	
Expiry date:		2021-06-30	
Coagulometers		MASTER CURVE	
		µg(FEU)/mL	ΔOD
Optical	Yumizen G400 DDi/G200	4,00	0,1058
		3,00	0,0792
		2,00	0,0525
		1,00	0,0262
		0,50	0,0133
		0,25	0,0070
		Slope value	N/A
Important notice			
<p>For Yumizen G400 DDi and Yumizen G200, make sure that the end time is set to 180 seconds in the D-Dimer setup and the curve type Lin-Lin p-p is selected in Calibration > D-Dimer > ug FEU/ml. For more information, please contact your local HORIBA Medical representative.</p>			