



Yumizen G PT 5 for Yumizen G Line



A 4 9 0 0 6 0 9 3 6 3 3 8 0 6 2 2 A



A 5 1 0 0 1 2 6 5 0 0 1 6 7 3 3 0 A



A 6 2 0 9 2 5 0 2 4 9 1 2 6 1 1 9 A



A 7 1 6 7 A

Reference No.:		1300036338	
Lot No.:		900609 (rev. 1)	
Expiry date:		2022-06-30	
Coagulometers		ISI	MASTER CURVE
Optical	Yumizen G Line		%
		100	12,6
		50	16,7
		33	20,9
		25	24,9
		MNPT	12,6
Important notice			
<ul style="list-style-type: none"> - The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT. - The efficiency of the reconstitution of the reagent may be approved by performing it at 37°C. 			