



Yumizen G DDi 2 for YUMIZEN G1500/G1550



D 4900912363910921D



D 5400263030020402D



D 6001265100055505D



D 7002800250155166D

Reference No.:		1300036391	
Lot No.:		900912 (rev. 1)	
Expiry date:		2021-09-30	
Coagulometers		MASTER CURVE	
Optical	Yumizen G1500/G1550	µg(FEU)/mL	ΔOD
		4,00	0,2630
		3,00	0,2040
		2,00	0,1265
		1,00	0,0555
		0,50	0,0280
		0,25	0,0155
Slope value		0,07	

Yumizen G DDi 2 for YUMIZEN G800



D 4900912363910921D



D 5400252530019452D



D 6001190100051005D



D 7002450250130164D

Reference No.:		1300036391	
Lot No.:		900912 (rev. 1)	
Expiry date:		2021-09-30	
Coagulometers		MASTER CURVE	
Optical	Yumizen G800	µg(FEU)/mL	ΔOD
		4,00	0,2525
		3,00	0,1945
		2,00	0,1190
		1,00	0,0510
		0,50	0,0245
		0,25	0,0130
Slope value		0,07	



Yumizen G DDi 2 for YUMIZEN G400 DDi/G200

Reference No.:		1300036391	
Lot No.:		900912 (rev. 1)	
Expiry date:		2021-09-30	
Coagulometers		MASTER CURVE	
		µg(FEU)/mL	ΔOD
Optical	Yumizen G400 DDi/G200	4,00	0,1160
		3,00	0,0880
		2,00	0,0570
		1,00	0,0270
		0,50	0,0130
		0,25	0,0065
		Slope value	N/A
Important notice			
<p>For Yumizen G400 DDi and Yumizen G200, make sure that the end time is set to 180 seconds in the D-Dimer setup and the curve type Lin-Lin p-p is selected in Calibration > D-Dimer > ug FEU/ml. For more information, please contact your local HORIBA Medical representative.</p>			