



LAtrol™
Normal Control Plasma
for YUMIZEN G800

Reference No.:		1300079435
Lot No.:		200626
Expiry date:		26/06/2022
Optical	Yumizen G800	
Parameter	Unit	Range
DRVVT	sec	42.4 – 54.4
DRVVC	sec	37.1 – 48.9
Normalized Ratio		<1.16

LAtrol™
Normal Control Plasma
for YUMIZEN
G1500/G1550

Reference No.:		1300079435
Lot No.:		200626
Expiry date:		26/06/2022
Optical	Yumizen G1500/G1550	
Parameter	Unit	Range
DRVVT	sec	ND
DRVVC	sec	ND
Normalized Ratio		<1.16

Important notice

Each laboratory should establish its own DRVVT and DRVVC normal reference range (mean value ± 2 S.D.).

The Normalized Ratio is calculated as follows:

$$\text{Normalized Ratio} = \frac{\text{DRVVT (sample)} / \text{DRVVT (mean clot time test)}}{\text{DRVVC (sample)} / \text{DRVVC (mean clot time confirm)}}$$