

**Yumizen G CTRL DDi I  
 for YUMIZEN  
 G1500/G1550**



D 4 9 0 0 5 1 6 4 6 4 1 4 0 5 2 2 D



D 5 0 1 9 0 6 9 0 8 3 D

Reference No.:		1300036414
Lot No.:		900516 (rev. 2)
Expiry date:		2022-05-31
Optical	Yumizen G1500/G1550	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	0,19 – 0,69

**Yumizen G CTRL DDi I  
 for YUMIZEN G800**



D 4 9 0 0 5 1 6 4 6 4 1 4 0 5 2 2 D



D 5 0 2 0 0 7 0 0 6 7 D

Reference No.:		1300036414
Lot No.:		900516 (rev. 2)
Expiry date:		2022-05-31
Optical	Yumizen G800	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	0,20 – 0,70

**Yumizen G CTRL DDi I  
 for YUMIZEN  
 G400 DDi/G200**

Reference No.:		1300036414
Lot No.:		900516 (rev. 2)
Expiry date:		2022-05-31
Optical	Yumizen G400 DDi/G200	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	0,13 – 0,63

**Yumizen G CTRL DDi II  
 for YUMIZEN  
 G1500/G1550**



Reference No.:		1300036414
Lot No.:		900516 (rev. 2)
Expiry date:		2022-05-31
Optical	Yumizen G1500/G1550	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	1,52 – 2,54

**Yumizen G CTRL DDi II  
 for YUMIZEN G800**



Reference No.:		1300036414
Lot No.:		900516 (rev. 2)
Expiry date:		2022-05-31
Optical	Yumizen G800	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	1,58 – 2,63

**Yumizen G CTRL DDi II  
 for YUMIZEN  
 G400 DDi/G200**

Reference No.:		1300036414
Lot No.:		900516 (rev. 2)
Expiry date:		2022-05-31
Optical	Yumizen G400 DDi/G200	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	1,43 – 2,43