



Yumizen G DDi 2 for YUMIZEN G1500/G1550



D 4901106363911121D



D 5400263930020852D



D 6001340100059005D



D 7002900250155167D

Reference No.:		1300036391	
Lot No.:		901106 (rev. 1)	
Expiry date:		2021-11-30	
Coagulometers		MASTER CURVE	
Optical	Yumizen G1500/G1550	µg(FEU)/mL	ΔOD
		4,00	0,2639
		3,00	0,2085
		2,00	0,1340
		1,00	0,0590
		0,50	0,0290
		0,25	0,0155
Slope value		0,07	

Yumizen G DDi 2 for YUMIZEN G800



D 4901106363911121D



D 5400243530018602D



D 6001190100054505D



D 7002650250140160D

Reference No.:		1300036391	
Lot No.:		901106 (rev. 1)	
Expiry date:		2021-11-30	
Coagulometers		MASTER CURVE	
Optical	Yumizen G800	µg(FEU)/mL	ΔOD
		4,00	0,2435
		3,00	0,1860
		2,00	0,1190
		1,00	0,0545
		0,50	0,0265
		0,25	0,0140
Slope value		0,07	



Yumizen G DDi 2 for YUMIZEN G400 DDi/G200

Reference No.:		1300036391	
Lot No.:		901106 (rev. 1)	
Expiry date:		2021-11-30	
Coagulometers		MASTER CURVE	
		µg(FEU)/mL	ΔOD
Optical	Yumizen G400 DDi/G200	4,00	0,1064
		3,00	0,0840
		2,00	0,0553
		1,00	0,0254
		0,50	0,0125
		0,25	0,0068
		Slope value	0,04
Important notice			
<p>For Yumizen G400 DDi and Yumizen G200, make sure that the end time is set to 180 seconds in the D-Dimer setup and the curve type Lin-Lin p-p is selected in Calibration > D-Dimer > ug FEU/ml. For more information, please contact your local HORIBA Medical representative.</p>			