



Yumizen G FIB 5 for Yumizen G1500/G1550 and Yumizen G800



D 4 9 0 0 4 0 4 3 6 3 8 4 0 4 2 2 D



D 5 4 0 0 0 6 4 3 0 0 0 8 0 2 0 0 D



D 6 1 0 9 1 0 0 1 8 5 0 0 0 0 0 0 D



D 7 1 2 3 D

Reference No.:		1300036384	
Lot No.:		900404 (rev. 2)	
Expiry date:		2022-04-30	
Coagulometers		MASTER CURVE	
Optical	Yumizen G1500/G1550/G800	g/L	sec
		4	6,4
		3	8,0
		2	10,9
		1	18,5
Dilution		10x	

Yumizen G FIB 5 for Yumizen G400 DDi/G200



D 4 9 0 0 4 0 4 3 6 3 8 4 0 4 2 2 D



D 5 4 0 0 0 6 4 3 0 0 0 8 1 2 0 0 D



D 6 1 1 2 1 0 0 1 9 5 0 0 0 0 0 0 D



D 7 1 1 9 D

Reference No.:		1300036384	
Lot No.:		900404 (rev. 2)	
Expiry date:		2022-04-30	
Coagulometers		MASTER CURVE	
Optical	Yumizen G400 DDi/G200	g/L	sec
		4	6,4
		3	8,1
		2	11,2
		1	19,5
Dilution		10x	