



# Yumizen G FIB 5 for Yumizen G1500/G1550 and Yumizen G800



D 4900808363840822D



D 5400070300086200D



D 6114100186000000D



D 7135D

<b>Reference No.:</b>		<b>1300036384</b>	
<b>Lot No.:</b>		<b>900808 (rev. 2)</b>	
<b>Expiry date:</b>		<b>2022-08-31</b>	
<b>Coagulometers</b>		<b>MASTER CURVE</b>	
<b>Optical</b>	Yumizen G1500/G1550/G800	<b>g/L</b>	<b>sec</b>
		4	7,0
		3	8,6
		2	11,4
		1	18,6
<b>Dilution</b>		<b>10x</b>	

# Yumizen G FIB 5 for Yumizen G400 DDi/G200



D 4900808363840822D



D 5400064300082200D



D 6115100204000000D



D 7126D

<b>Reference No.:</b>		<b>1300036384</b>	
<b>Lot No.:</b>		<b>900808 (rev. 2)</b>	
<b>Expiry date:</b>		<b>2022-08-31</b>	
<b>Coagulometers</b>		<b>MASTER CURVE</b>	
<b>Optical</b>	Yumizen G400 DDi/G200	<b>g/L</b>	<b>sec</b>
		4	6,4
		3	8,2
		2	11,5
		1	20,4
<b>Dilution</b>		<b>10x</b>	