



# Yumizen G FIB 2 for Yumizen G Line



|                       |                       |                        |             |
|-----------------------|-----------------------|------------------------|-------------|
| <b>Reference No.:</b> |                       | <b>1300036383</b>      |             |
| <b>Lot No.:</b>       |                       | <b>901213 (rev. 1)</b> |             |
| <b>Expiry date:</b>   |                       | <b>2022-12-31</b>      |             |
| <b>Coagulometers</b>  |                       | <b>MASTER CURVE</b>    |             |
| <b>Optical</b>        | <b>Yumizen G Line</b> | <b>g/L</b>             | <b>sec</b>  |
|                       |                       | <b>4</b>               | <b>7,0</b>  |
|                       |                       | <b>3</b>               | <b>8,6</b>  |
|                       |                       | <b>2</b>               | <b>11,6</b> |
|                       |                       | <b>1</b>               | <b>19,4</b> |
| <b>Dilution</b>       |                       | <b>10x</b>             |             |