



Yumizen G CTRL DDi I for YUMIZEN G1500/G1550



D 4 9 1 0 3 2 3 4 6 4 1 4 0 3 2 3 D



D 5 0 1 9 0 6 9 0 7 9 D

Reference No.:		1300036414
Lot No.:		910323 (rev. 1)
Expiry date:		2023-03-31
Optical	Yumizen G1500/G1550	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	0,19 – 0,69

Yumizen G CTRL DDi I for YUMIZEN G800



D 4 9 1 0 3 2 3 4 6 4 1 4 0 3 2 3 D



D 5 0 1 9 0 6 9 0 7 9 D

Reference No.:		1300036414
Lot No.:		910323 (rev. 1)
Expiry date:		2023-03-31
Optical	Yumizen G800	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	0,19 – 0,69

Yumizen G CTRL DDi I for YUMIZEN G400 DDi/G200

Reference No.:		1300036414
Lot No.:		910323 (rev. 1)
Expiry date:		2023-03-31
Optical	Yumizen G400 DDi/G200	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	0,16 – 0,66



**Yumizen G CTRL DDi II
for YUMIZEN
G1500/G1550**



D 4 9 1 0 3 2 3 5 6 4 1 4 0 3 2 3 D



D 5 1 3 1 2 3 7 0 7 2 D

Reference No.:		1300036414
Lot No.:		910323 (rev. 1)
Expiry date:		2023-03-31
Optical	Yumizen G1500/G1550	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	1,31 – 2,37

**Yumizen G CTRL DDi II
for YUMIZEN G800**



D 4 9 1 0 3 2 3 5 6 4 1 4 0 3 2 3 D



D 5 1 3 3 2 4 0 0 6 8 D

Reference No.:		1300036414
Lot No.:		910323 (rev. 1)
Expiry date:		2023-03-31
Optical	Yumizen G800	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	1,33 – 2,40

**Yumizen G CTRL DDi II
for YUMIZEN
G400 DDi/G200**

Reference No.:		1300036414
Lot No.:		910323 (rev. 1)
Expiry date:		2023-03-31
Optical	Yumizen G400 DDi/G200	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	1,37 – 2,47