



**Yumizen G CTRL DDi I
 for YUMIZEN
 G1500/G1550**



D 4 9 1 0 4 2 9 4 6 4 1 4 0 4 2 3 D



D 5 0 2 0 0 7 0 0 7 1 D

Reference No.:		1300036414
Lot No.:		910429 (rev. 1)
Expiry date:		2023-04-30
Optical	Yumizen G1500/G1550	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	0,20 – 0,70

**Yumizen G CTRL DDi I
 for YUMIZEN G800**



D 4 9 1 0 4 2 9 4 6 4 1 4 0 4 2 3 D



D 5 0 2 0 0 7 0 0 7 1 D

Reference No.:		1300036414
Lot No.:		910429 (rev. 1)
Expiry date:		2023-04-30
Optical	Yumizen G800	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	0,20 – 0,70

**Yumizen G CTRL DDi I
 for YUMIZEN
 G400 DDi/G200**

Reference No.:		1300036414
Lot No.:		910429 (rev. 1)
Expiry date:		2023-04-30
Optical	Yumizen G400 DDi/G200	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	0,13 – 0,63



**Yumizen G CTRL DDi II
for YUMIZEN
G1500/G1550**



D 4 9 1 0 4 2 9 5 6 4 1 4 0 4 2 3 D



D 5 1 3 6 2 4 5 0 8 4 D

Reference No.:		1300036414
Lot No.:		910429 (rev. 1)
Expiry date:		2023-04-30
Optical	Yumizen G1500/G1550	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	1,36 – 2,45

**Yumizen G CTRL DDi II
for YUMIZEN G800**



D 4 9 1 0 4 2 9 5 6 4 1 4 0 4 2 3 D



D 5 1 3 8 2 4 9 0 9 0 D

Reference No.:		1300036414
Lot No.:		910429 (rev. 1)
Expiry date:		2023-04-30
Optical	Yumizen G800	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	1,38 – 2,49

**Yumizen G CTRL DDi II
for YUMIZEN
G400 DDi/G200**

Reference No.:		1300036414
Lot No.:		910429 (rev. 1)
Expiry date:		2023-04-30
Optical	Yumizen G400 DDi/G200	
Parameter	Unit	Range
D-Dimer	µg(FEU)/mL	1,27 – 2,28