



# Yumizen G PT Reco 10 for Yumizen G Line



A 4 9 9 1 2 0 3 3 6 3 7 6 1 2 2 1 A



A 5 1 0 0 1 0 9 5 0 0 1 5 7 3 3 0 A



A 6 2 0 6 2 5 0 2 5 3 1 0 9 1 0 0 A



A 7 1 4 8 A

Reference No.:		1300036376		
Lot No.:		991203		
Expiry date:		2021-12-31		
Coagulometers		ISI	MASTER CURVE	
Optical	Yumizen G Line	1,00	%	SEC
			100	10,9
			50	15,7
			33	20,6
			25	25,3
MNPT				10,9
<b>Important notice</b>				
<p>The Mean Normal Prothrombin Time (MNPT) depends on the population, race, gender, sampling tube, etc. Our value, that is identical with the 100% point of the calibration curve is for information only. According to the CLSI every laboratory should determine its own MNPT.</p> <p>For Yumizen G200, Yumizen G400 and Yumizen G400 DDi, the <b>Min. Time</b> must be set to 7 seconds and the <b>Lag time</b> to 5 seconds in the PT setup. Please refer to the Yumizen G PT Reco 10 reagent notice.</p>				