



# Yumizen G FIB 2 for Yumizen G1500/G1550 and Yumizen G800



Reference No.:		1300036383	
Lot No.:		910418 (rev. 1)	
Expiry date:		2023-04-30	
Coagulometers		MASTER CURVE	
Optical	Yumizen G1500/G1550/G800	g/L	sec
		4	8,1
		3	9,8
		2	12,9
		1	20,4
		Dilution	10x

# Yumizen G FIB 2 for Yumizen G400 DDi/G200



Reference No.:		1300036383	
Lot No.:		910418 (rev. 1)	
Expiry date:		2023-04-30	
Coagulometers		MASTER CURVE	
Optical	Yumizen G400 DDi/G200	g/L	sec
		4	8,8
		3	10,9
		2	14,5
		1	23,8
		Dilution	10x