

# Yumizen G FIB 2 for Yumizen G1500/G1550 and Yumizen G800



Reference No.:		1300036383	
Lot No.:		910604 (rev. 1)	
Expiry date:		2023-06-30	
Coagulometers		MASTER CURVE	
Optical	Yumizen G1500/G1550/G800	g/L	sec
		4	7,8
		3	9,5
		2	12,3
		1	19,5
		Dilution	10x

# Yumizen G FIB 2 for Yumizen G400 DDi/G200



Reference No.:		1300036383	
Lot No.:		910604 (rev. 1)	
Expiry date:		2023-06-30	
Coagulometers		MASTER CURVE	
Optical	Yumizen G400 DDi/G200	g/L	sec
		4	8,2
		3	10,2
		2	13,8
		1	23,0
		Dilution	10x