



Yumizen G FIB 5 for Yumizen G1500/G1550 and Yumizen G800



| | | | |
|-----------------------|--------------------------|------------------------|-------------|
| Reference No.: | | 1300036384 | |
| Lot No.: | | 910304 (rev. 1) | |
| Expiry date: | | 2023-03-31 | |
| Coagulometers | | MASTER CURVE | |
| | | g/L | sec |
| Optical | Yumizen G1500/G1550/G800 | 4 | 8,0 |
| | | 3 | 9,9 |
| | | 2 | 13,3 |
| | | 1 | 22,0 |
| | | Dilution | 10x |

Yumizen G FIB 5 for Yumizen G400 DDi/G200



| | | | |
|-----------------------|-----------------------|------------------------|-------------|
| Reference No.: | | 1300036384 | |
| Lot No.: | | 910304 (rev. 1) | |
| Expiry date: | | 2023-03-31 | |
| Coagulometers | | MASTER CURVE | |
| | | g/L | sec |
| Optical | Yumizen G400 DDi/G200 | 4 | 8,8 |
| | | 3 | 10,9 |
| | | 2 | 14,8 |
| | | 1 | 25,1 |
| | | Dilution | 10x |