

Yumizen G FIB 2 for Yumizen G1500/G1550 and Yumizen G800



Reference No.:		1300036383	
Lot No.:		910707 (rev. 1)	
Expiry date:		2023-07-31	
Coagulometers		MASTER CURVE	
Optical	Yumizen G1500/G1550/G800	g/L	sec
		4	8,2
		3	9,8
		2	12,6
		1	19,2
		Dilution	10x

Yumizen G FIB 2 for Yumizen G400 DDi/G200



Reference No.:		1300036383	
Lot No.:		910707 (rev. 1)	
Expiry date:		2023-07-31	
Coagulometers		MASTER CURVE	
Optical	Yumizen G400 DDi/G200	g/L	sec
		4	8,6
		3	10,5
		2	13,9
		1	22,3
		Dilution	10x