



# Yumizen G FIB 5 for Yumizen G1500/G1550 and Yumizen G800



D 4 9 1 0 7 0 8 3 6 3 8 4 0 7 2 3 D



D 5 4 0 0 0 8 3 3 0 0 0 9 9 2 0 0 D



D 6 1 2 7 1 0 0 1 9 5 0 0 0 0 0 0 D



D 7 1 4 7 D

<b>Reference No.:</b>		<b>1300036384</b>	
<b>Lot No.:</b>		<b>910708 (rev. 1)</b>	
<b>Expiry date:</b>		<b>2023-07-31</b>	
<b>Coagulometers</b>		<b>MASTER CURVE</b>	
<b>Optical</b>	<b>Yumizen G1500/G1550/G800</b>	<b>g/L</b>	<b>sec</b>
		<b>4</b>	<b>8,3</b>
		<b>3</b>	<b>9,9</b>
		<b>2</b>	<b>12,7</b>
		<b>1</b>	<b>19,5</b>
	<b>Dilution</b>	<b>10x</b>	

# Yumizen G FIB 5 for Yumizen G400 DDi/G200



D 4 9 1 0 7 0 8 3 6 3 8 4 0 7 2 3 D



D 5 4 0 0 0 8 6 3 0 0 1 0 5 2 0 0 D



D 6 1 3 8 1 0 0 2 2 3 0 0 0 0 0 0 D



D 7 1 3 2 D

<b>Reference No.:</b>		<b>1300036384</b>	
<b>Lot No.:</b>		<b>910708 (rev. 1)</b>	
<b>Expiry date:</b>		<b>2023-07-31</b>	
<b>Coagulometers</b>		<b>MASTER CURVE</b>	
<b>Optical</b>	<b>Yumizen G400 DDi/G200</b>	<b>g/L</b>	<b>sec</b>
		<b>4</b>	<b>8,6</b>
		<b>3</b>	<b>10,5</b>
		<b>2</b>	<b>13,8</b>
		<b>1</b>	<b>22,3</b>
	<b>Dilution</b>	<b>10x</b>	