

LAtrol™
Normal Control Plasma
for YUMIZEN G800

Reference No.:		1300079435
Lot No.:		210811_Rev2
Expiry date:		08/11/2023
Optical	Yumizen G800	
Parameter	Unit	Range
DRVVT	sec	38.5 – 51.9
DRVVC	sec	27.6 – 37.1
Normalized Ratio		<1.16

LAtrol™
Normal Control Plasma
for YUMIZEN
G1500/G1550

Reference No.:		1300079435
Lot No.:		210811
Expiry date:		08/11/2023
Optical	Yumizen G1500/G1550	
Parameter	Unit	Range
DRVVT	sec	-
DRVVC	sec	-
Normalized Ratio		<1.16

Important notice

Each laboratory should establish its own DRVVT and DRVVC normal reference range (mean value ± 2 S.D.).

The Normalized Ratio is calculated as follows:

$$\text{Normalized Ratio} = \frac{\text{DRVVT (sample)} / \text{DRVVT (mean clot time test)}}{\text{DRVVC (sample)} / \text{DRVVC (mean clot time confirm)}}$$