



# Yumizen G DDi 2 for YUMIZEN G1500/G1550



D 4911129363911122D



D 5400235030017602D



D 6001105100050505D



D 7002450250125155D

<b>Reference No.:</b>		<b>1300036391</b>	
<b>Lot No.:</b>		<b>911129 (rev. 1)</b>	
<b>Expiry date:</b>		<b>2022-11-30</b>	
<b>Coagulometers</b>		<b>MASTER CURVE</b>	
<b>Optical</b>	<b>Yumizen G1500/G1550</b>	<b>µg(FEU)/mL</b>	<b>ΔOD</b>
		4,00	0,2350
		3,00	0,1760
		2,00	0,1105
		1,00	0,0505
		0,50	0,0245
		0,25	0,0125
<b>Slope value</b>		<b>0,05</b>	

# Yumizen G DDi 2 for YUMIZEN G800



D 4911129363911122D



D 5400222530016702D



D 6001030100046005D



D 7002160250110145D

<b>Reference No.:</b>		<b>1300036391</b>	
<b>Lot No.:</b>		<b>911129 (rev. 1)</b>	
<b>Expiry date:</b>		<b>2022-11-30</b>	
<b>Coagulometers</b>		<b>MASTER CURVE</b>	
<b>Optical</b>	<b>Yumizen G800</b>	<b>µg(FEU)/mL</b>	<b>ΔOD</b>
		4,00	0,2225
		3,00	0,1670
		2,00	0,1030
		1,00	0,0460
		0,50	0,0216
		0,25	0,0110
<b>Slope value</b>		<b>0,06</b>	



# ***Yumizen G DDi 2*** ***for YUMIZEN*** ***G400 DDi/G200***

<b>Reference No.:</b>		<b>1300036391</b>	
<b>Lot No.:</b>		<b>911129 (rev. 1)</b>	
<b>Expiry date:</b>		<b>2022-11-30</b>	
<b>Coagulometers</b>		<b>MASTER CURVE</b>	
		<b>µg(FEU)/mL</b>	<b>ΔOD</b>
<b>Optical</b>	<b>Yumizen G400 DDi/G200</b>	<b>4,00</b>	<b>0,0898</b>
		<b>3,00</b>	<b>0,0665</b>
		<b>2,00</b>	<b>0,0428</b>
		<b>1,00</b>	<b>0,0203</b>
		<b>0,50</b>	<b>0,0100</b>
		<b>0,25</b>	<b>0,0050</b>
		<b>Slope value</b>	<b>0,02</b>
<b>Important notice</b>			
<p>For Yumizen G400 DDi and Yumizen G200, make sure that the end time is set to 180 seconds in the D-Dimer setup and the curve type <b>Lin-Lin p-p</b> is selected in <b>Calibration &gt; D-Dimer &gt; ug FEU/ml</b>.          For more information, please contact your local HORIBA Medical representative.</p>			