



**HORIBA ABX SAS**  
 Parc Euromédecine  
 Rue du Caducée  
 B.P. 7290  
 34184 Montpellier cedex 4  
 FRANCE  
 www.horiba.com



A03A00013C

# Yumizen G FIB 5 for Yumizen G1500/G1550 and Yumizen G800



D 4 9 2 0 1 1 0 3 6 3 8 4 0 1 2 4 D



D 5 4 0 0 0 8 2 3 0 0 0 9 7 2 0 0 D



D 6 1 2 5 1 0 0 1 9 0 0 0 0 0 0 0 D



D 7 1 2 0 D

<b>Reference No.:</b>		<b>1300036384</b>	
<b>Lot No.:</b>		<b>920110 (rev. 1)</b>	
<b>Expiry date:</b>		<b>2024-01-31</b>	
<b>Coagulometers</b>		<b>MASTER CURVE</b>	
<b>Optical</b>	Yumizen G1500/G1550/G800	<b>g/L</b>	<b>sec</b>
		4	8,2
		3	9,7
		2	12,5
		1	19,0
	<b>Dilution</b>	<b>10x</b>	

# Yumizen G FIB 5 for Yumizen G400 DDi/G200



D 4 9 2 0 1 1 0 3 6 3 8 4 0 1 2 4 D



D 5 4 0 0 0 8 6 3 0 0 1 0 5 2 0 0 D



D 6 1 4 1 1 0 0 2 3 3 0 0 0 0 0 0 D



D 7 1 1 0 D

<b>Reference No.:</b>		<b>1300036384</b>	
<b>Lot No.:</b>		<b>920110 (rev. 1)</b>	
<b>Expiry date:</b>		<b>2024-01-31</b>	
<b>Coagulometers</b>		<b>MASTER CURVE</b>	
<b>Optical</b>	Yumizen G400 DDi/G200	<b>g/L</b>	<b>sec</b>
		4	8,6
		3	10,5
		2	14,1
		1	23,3
	<b>Dilution</b>	<b>10x</b>	