

ABX Minocal



(Exp.) 2016-10-05

(YYYY - MM - DD)

LOT		CX 402		CAL								
		Rev 1				ABX Lyse						
PARAMETRES PARAMETERS	UNITES UNITS	MICROS		MICROS		MICROS		PENTRA		PENTRA		TOLERANCES TOLERANCE
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P 60 18P ADVIA 60						
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.4	9.3	9.3	9.3	9.4						± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52	4.48	4.53	4.51	4.50						± 0.06
HB HGB	g/dl	13.7	13.7	13.7	13.6	13.6						± 0.2
	g/l	137	137	137	136	136						± 2
	mmol/l	8.51	8.51	8.51	8.45	8.45						± 0.12
HT HCT	%	38.8	38.5	39.0	38.8	38.6						± 1.0
	l/l	0.388	0.385	0.390	0.388	0.386						± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	243	240	242	245	241						± 10
VPM MPV	µm ³ ; fl	9.1	9.0	N/A	9.2	8.2						± 0.5

Ref: TEMP-0387 Rev.42 FRONT / RECTO 1219930090 -A

LOT		CX 402		CAL		 (Exp.) 2016-10-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		MICROS	MICROS	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	PENTRA	
		EMI CRP		MS CRP	XLR	DX NEXUS DF NEXUS	MS60	120 - 120 RETIC DX120 - DF120	60 60 C+ ES60	80 XL80	
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.4		9.4	9.3	9.5	9.3	9.5	9.3	9.2	± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.48		4.48	4.53	4.53	4.46	4.53	4.52	4.49	± 0.06
HB HGB	g/dl	13.5		13.4	13.6	13.5	13.5	13.5	13.6	13.6	± 0.2
	g/l	135		134	136	135	135	135	136	136	± 2
	mmol/l	8.38		8.32	8.45	8.38	8.38	8.38	8.45	8.45	± 0.12
HT HCT	%	38.1		36.5	37.1	37.4	36.1	37.4	37.1	37.0	± 1.0
	l/l	0.381		0.365	0.371	0.374	0.361	0.374	0.371	0.370	± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	246		230	240	245	240	245	243	240	± 10
VPM MPV	µm ³ ; fl	9.0		N/A	N/A	N/A	N/A	N/A	N/A	N/A	± 0.5

Ref: TEMP-0387 Rev.42 BACK / VERSO 1219930090 -A